

Post Mortem Tissue Distribution of Synthetic Cathinones

Lindsay Glicksberg, PhD^{1,2}; Ruth Winecker, PhD³; Caitlin Miller, MS⁴;
Sarah Kerrigan, PhD¹

¹Department of Forensic Sciences, Sam Houston State University;

²Dallas County Southwestern Institute of Forensic Sciences;

³North Carolina Department of Health and Human Services;

⁴LA County Department of Medical Examiner-Coroner



FSF Emerging Forensic Scientist Award
Paper Presentation



Disclosure

- There is no real or apparent conflicts of interest related to the content of this presentation
- Products used:
 - Agilent Technologies 6530 Accurate-Mass Q-TOF LC/MS
- The authors declare no competing interest
- Project funded by NIH Award # 2013-R2-CX-K006



Postmortem Toxicology

Pre-analytical Factors

- Environment
- Postmortem Interval
- Drug Properties

Further Complications

- Drug Instability
- Postmortem Redistribution (PMR)

Challenging Samples

- Sample Quality
- Sample Collection
- Alternative Samples

Interpretation Challenges

- Changes in concentration
- Site to site differences

Factors Influencing PMR

- Drug Properties
 - Volume of Distribution ($>3-4$)
 - pK_a (>7)
 - Lipophilicity
 - Physio-Chemical Changes
 - Cell death
 - Blood acidification
 - Drug Release
 - Organs
 - Protein bound
 - Blood Movement
 - Putrefactive processes
 - Physical
 - Bacteria
 - External
 - GI
 - Environment
 - Postmortem Interval
 - Incomplete Distribution
-
- Pre-analytical factors
 - Case, time, and site dependent variables
 - Basic and lipophilic drugs more susceptible to redistribution




Assessing PMR

- Antemortem to postmortem specimen comparison
 - Possible in hospital deaths
- Cardiac to Peripheral blood (C/P) ratio
 - $C/P > 1$
- Liver to Peripheral blood (L/P) ratio
 - $L/P > 20$
 - $L/P < 5$
- More reliable with large population of cases
 - Case to case variability


General Structure of Cathinone

- Arylamino ketones
- R groups represent positions that can be substituted to create various synthetic cathinones




General Structure of Cathinone

- Arylaminoketones
- R groups represent positions that can be substituted to create various synthetic cathinones



Methamphetamine



Methcathinone

Synthetic Cathinone PMR

C/P—greater than 1; L/P—less than 5

- Individual cases
- Small sample population

Cathinone	C/P	L/P	References
Butylone	-	1.7	Rojek, 2012
Ethylone	1.0	3.6	McIntyre, 2014
MDPV	0.7 – 1.7 (1.3, n=6)	2.2 – 23 (8.9, n=7)	Marinetti, 2013; Wyman, 2013; Kesha, 2013
Methylone	1.0 – 2.1 (1.3, n=5)	1.6 – 3.2 (2.6, n=4)	Shimomura, 2016; Pearson, 2012; McIntyre, 2013; Cawrse, 2012; Barrios, 2016
Pentedrone	-	1.1	Sykutera, 2015
α -PVP	1.45	1.0 – 2.9 (1.7, n=3)	Potocka, 2017; Hasegawa, 2014; Sykutera, 2015
Pyrovalerone	1.4	3.0	Marinetti, 2013



Statement of Problem

Synthetic cathinones associated with fatal intoxications

Postmortem results require careful interpretation

Goal: Assess tissue distribution and postmortem redistribution of synthetic cathinones

Participating Laboratories

North Carolina Department of Health and Human Services, Office of the Chief Medical Examiner

- ethylone (n=8)
- methylone (n=14)
- α -PVP (n=16)
- pentedrone (n=1)
- 4-MEC (n=1)
- pentylone (n=1)
- methedrone (n=1)
- buytlone (n=1)
- MDPV (n=4)

LA County Department of Medical Examiner-Coroner

- Ethylone (n=3)

50 cases, 139 specimens total



Specimens

Central blood sources:

- Vena cava
- Right chest cavity
- Heart blood
- Central blood
- Aorta blood

Vitreous humor

Stomach contents

Urine

Liver

Peripheral blood sources:

- Subclavian vessel
- Peripheral blood
- Jugular vein
- Iliac vein
- Femoral vessel

Optimized Extraction

Bovine blood (2 mL) (25 ng/mL I.S.)

Protein precipitation (4 mL ACN)

Centrifuge (5 min, 4,000 RPM)

1.0 M PO₄ buffer, pH 6 (6 mL)

Drug free urine (1 mL)
(25 ng/mL I.S.)

1.0 M PO₄ buffer, pH 6 (2 mL)

Add to SPE cartridges

Aqueous Washes

- DI water (1 mL)
- 1 M acetic acid (1 mL)

Vacuum Dry (5 min)

Organic Washes

- Hexane (1 mL)
- Ethyl acetate (1 mL)
- Methanol (1 mL)
- Dichloromethane (1 mL) (blood only)

Elute: 2% NH₄OH in 95 DCM: 5 IPA (1 mL)

Evaporation (50°C, N₂)

Reconstitute 50:50 MP A:B (25 µL)

Specimen Extraction

Liver

- Homogenization
 - Bead Ruptor 12
 - One part tissue (0.5 g) to two parts deionized water
 - 0.5 mL homogenate fortified and extracted
- Blood protocol
 - 2 mL ACN protein precipitation
 - 3 mL phosphate buffer

Other Matrices

- Vitreous, stomach contents
 - Urine protocol



LC-Q/TOF MS Conditions

Agilent Technologies 6530 Accurate-Mass Q-TOF LC/MS

LC Separation

- Poroshell 120 EC-C18 Column (2.1x100 mm, 2.7 μ m particle size)
- Mobile Phase A: 0.1% FA in diH₂O
- Mobile Phase B: 0.1% FA in ACN
- Flow Rate: 0.40 mL/min
- LC Gradient:
 - 96% A to 5 min, 90% A until 11 min, 60% A for 1 min, 0% A to equilibrate the column

Q/TOF Parameters

- Gas Temperature: 200°C
- Gas Flow Rate: 13 L/min
- Sheath Gas Temperature: 250°C
- Sheath Gas Flow Rate: 12 L/min
- Nebulizer Pressure: 20 psig

Mass Spectrometry

- Capillary Voltage: 4000 V
- Fragmentor Voltage: 150 V
- Nozzle Voltage: 0 V
- Collision Energy: 30 eV, 20 eV
- MS Scan Rate: 8 spectra/sec
- MS/MS Scan Rate: 3 spectra/sec
- MS Scan Range: 40-1000 m/z
- ESI Mode: Positive

Acquisition

- Minimum of two ion transitions per drug
- Run Time: 13 minutes

Validation Summary

Evaluated according to Scientific Working Group of Forensic Toxicology (SWGTOX)
Standard Practices for Method Validation

Parameter	Blood	Urine
Calibration Model	Quadratic, 1/x	
Limit of Detection	1 – 5 ng/mL	0.25 – 5 ng/mL
Limit of Quantitation	1 – 5 ng/mL	0.25 – 5 ng/mL
Inter-assay Precision (n=15)	3 – 12%	2 – 12%
Intra-assay Precision (n=3)	0 – 17%	0 – 11%
Bias (n=15)	-7 – 11%	-3 – 12%
Accuracy (n=15)	93 – 111%	97 – 112%
Analytical Recovery (n=4)	81 – 91%	84 – 104%
Interference Studies	No qualitative interferences observed in either matrix	
Matrix Effects (n=10)	-15 – 3%	-21 – -1%
Dilution Integrity	2- and 4- fold	


Tissue Validation

Precision and Bias ($\pm 20\%$ CV)


- Low: 10 ng
- Mid: 25 ng
- High: 50 ng

Results

- All within 15%
- Exceptions:
 - 3-FMC (0.3-26.4% CV)
 - Ethcathinone (0.7-20.9% CV)
 - MDPBP (2.1-20.7% CV)
 - MPBP (0.7-24.7% CV)



Case #48
Ethylone
Positive



Cathinone-Positive Cases

50 positive cases

Cathinones identified:


- α -PVP (n=19)
- Methylone (n=18)
- Ethylone: (n=15)
- MDPV: (n=6)
- Pentylone: (n=3)
- Methedrone: (n=2)
- Butylone, 4-MEC, MDPBP: (n=1)

7 cases involved 2 or more cathinones


Vitreous (n=1), stomach contents (n=1)

Synthetic Cathinones Identified


Methedrone (n=2)




4-MEC (n=1)




Methylone (n=18)




Ethylone (n=15)




Pentylone (n=3)




Butylone (n=1)




α -PVP (n=19)



MDPV (n=6)



MDPBP (n=1)



Quantitative Results

Cathinone	Central Blood	Peripheral Blood	Liver	Urine
α -PVP	<2 – 1,090 ng/mL (n=17)	2 – 1,019 ng/mL (n=14)	<60 – 169 ng/g (n=5)	33 – 7,580 ng/mL (n=11)
Methylone	<2 – 202 ng/mL (n=19)	<2 – 28 ng/mL (n=7)	<60 – 1,347 ng/g (n=8)	2 – 38,064 ng/mL (n=11)
Ethylone	3 – 2,743 ng/mL (n=12)	<2 – 780 ng/mL (n=9)	<60 – 5,196 ng/g (n=6)	32 – 8,743 ng/mL (n=9)
MDPV	4 – 80 ng/mL (n=4)	3 – 80 ng/mL (n=2)	<60 – 223 ng/g (n=2)	4 – 5,210 ng/mL (n=3)
Pentylone	<5 – 323 ng/mL (n=2)	<5 – 160 ng/mL (n=2)		122 - >5,000 ng/mL (n=2)
Methedrone	79 ng/mL	70 ng/mL	720 ng/g	4 – 1,213 ng/mL (n=2)
Butylone	6 ng/mL	8 ng/mL	116 ng/g	934 ng/mL
4-MEC		57 ng/mL		
MDPBP				111 ng/mL

Postmortem Redistribution


Determine C/P and L/P ratios for 7 cathinones over 21 specimens

Cathinone	C/P	L/P
α -PVP	0.5 – 1.9 (1.1, n=9)	<1.4 (n=3)
MDPV	1.0 (n=1)	-
Methedrone	1.1 (n=1)	10 (n=1)
Butylone	0.7 (n=1)	14 (n=1)
Ethylone	0.5 – 9.2 (2.9, n=6)	0.2 – 20 (7.2, n=4)
Methylone	1.5 – 6.1 (4.0, n=5)	3.1 – 40 (n=2)
Pentylone	2.0 (n=1)	-

C/P Ratios

Determine C/P and L/P ratios for 7 cathinones

Cathinone	C/P	L/P
α -PVP	0.5 – 1.9 (1.1, n=9)	<1.4 (n=3)
MDPV	1.0 (n=1)	-
Methedrone	1.1 (n=1)	10 (n=1)
Butylone	0.7 (n=1)	14 (n=1)
Ethylone	0.5 – 9.2 (2.9, n=6)	0.2 – 20 (7.2, n=4)
Methylone	1.5 – 6.1 (4.0, n=5)	3.1 – 40 (n=2)
Pentylone	2.0 (n=1)	-



- Welch Test
 - methylone to ethylone
 - α -PVP to methylone
 - α -PVP to ethylone

Comparison to Literature

			Previously Published Values	
Cathinone	C/P	L/P	C/P	L/P
α -PVP	0.5 – 1.9 (1.1, n=9)	<1.4 (n=3)	1.5 (n=1)	1.1 – 2.9 (1.5, n=3)
MDPV	1.0 (n=1)	-	0.7 – 1.7 (1.3, n=6)	2.2 – 23 (8.9, n=7)
Methedrone	1.1 (n=1)	10 (n=1)	-	-
Butylone	0.7 (n=1)	14 (n=1)		1.7 (n=1)
Ethylone	0.5 – 9.2 (2.9, n=6)	0.2 – 20 (7.2, n=4)	1.0 (n=1)	3.6 (n=1)
Methylone	1.5 – 6.1 (4.0, n=5)	3.1 – 40 (n=2)	1.0 – 2.1 (1.3, n=5)	1.6 – 3.2 (2.6, n=4)
Pentylone	2.0 (n=1)	-	-	-

Comparison to Literature

			Previously Published Values		
Cathinone	C/P	L/P	C/P	L/P	Non Cathinones
α -PVP	0.5 – 1.9 (1.1, n=9)	<1.4 (n=3)	1.5 (n=1)	1.1 – 2.9 (1.5, n=3)	MDMA
MDPV	1.0 (n=1)	-	0.7 – 1.7 (1.3, n=6)	2.2 – 23 (8.9, n=7)	<u>C/P</u> : 0.9 – 4.6 (2.7, n=7)
Methedrone	1.1 (n=1)	10 (n=1)	-	-	<u>L/P</u> : 3.1 – 8.5 (6.5, n=5)
Butylone	0.7 (n=1)	14 (n=1)		1.7 (n=1)	Methamphetamine
Ethylone	0.5 – 9.2 (2.9, n=6)	0.2 – 20 (7.2, n=4)	1.0 (n=1)	3.6 (n=1)	<u>C/P</u> : 0.9 – 5.8 (2.3, n=39)
Methylone	1.5 – 6.1 (4.0, n=5)	3.1 – 40 (n=2)	1.0 – 2.1 (1.3, n=5)	1.6 – 3.2 (2.6, n=4)	<u>L/P</u> : 1.7 – 9.1 (5.5, n=19)
Pentylone	2.0 (n=1)	-	-	-	

Conclusions

- Cathinones may exhibit low to moderate PMR
- Highest C/P ratios observed for methylone and ethylone (secondary amine, methylenedioxy-substituted cathinones)
- Lowest C/P ratios observed for pyrrolidine-type cathinones

Limitations:

- Site and time dependent variables
- Further complicated by stability

Forensic Toxicology
<https://doi.org/10.1007/s11419-018-0403-3>



ORIGINAL ARTICLE



Postmortem distribution and redistribution of synthetic cathinones

Lindsay Glicksberg^{1,2} · Ruth Winecker³ · Caitlin Miller⁴ · Sarah Kerrigan¹

Received: 6 November 2017 / Accepted: 4 January 2018
© Japanese Association of Forensic Toxicology and Springer Japan KK, part of Springer Nature 2018

The header image features a dark blue background. On the right side, there is a large, detailed fingerprint. On the left side, there are faint, overlapping handwritten notes in a cursive script. The word 'Acknowledgements' is written in a large, bold, white sans-serif font across the center of the image.

Acknowledgements

This project was supported by Award No. 2013-R2-CX-K006 awarded by the National Institute of Justice, Office of Justice Programs, U.S. Department of Justice. The opinions, findings, and conclusions or recommendations expressed in this publication are those of the author(s) and do not necessarily reflect those of the Department of Justice.



FSF Emerging Forensic Scientist Award
Paper Presentation

Questions?

Lindsay Glicksberg
lcg012@shsu.edu

