Detection of Gasoline in Fire Debris Using Artificial Intelligence: Image Transformation of Gas Chromatography and Mass Spectrometry Data and Deep Learning

Ting-Yu Huang, MS*, and Jorn Chi-Chung Yu, PhD, ABC-CC


Department of Forensic Science Sam Houston State University Huntsville, TX 77340



Gasoline in arson fires

- Annual average of 52,260 arson fires were reported (2014-2018)
- Caused 400 civilian deaths, 950 civilian injuries, and \$815 million in direct property damage
- Gasoline is one of the most commonly-used ignitable liquids in arson fires
 - > Easy to obtain and transport
- Standard test method for identifying gasoline is ASTM E1618-19
 - Visual comparison of chromatograms
 - Extracted ion profiling
 - Target compound analysis








Deep learning (1/2)

- Convolutional neural network (CNN)
 - ➤ A neural network with multiple layers
 - Automatically extract features
 - Unstructured data (images)




Deep learning (2/2)

Advantages

- Capable of extract complex features
- Higher classification performance
- Predictive modeling
- Superior capability to classify images
- Wide applications in diagnosing diseases





Disadvantages

- Requires large-scale data collection
- High computational cost
- Transfer learning


 Can deep learning (transfer learning) be applied in gasoline detection if GC/MS data are presented in images?

Sample collection and preparation (1/3)

• Five brands of gasoline were collected in Huntsville, Texas (Brand A, B, C, D, and E)













Sample collection and preparation (2/3)

Serial dilution of the gasoline samples (n=8 dilution series in triplicates)






	Analyte concentration
N1	100 μg gasoline/20-mL HS vial
N2	50 μg gasoline/20-mL HS vial
N3	25 μg gasoline/20-mL HS vial
N4	12.5 μg gasoline/20-mL HS vial
N5	6.3 μg gasoline/20-mL HS vial
N6	3.1 μg gasoline/20-mL HS vial
N7	1.6 μg gasoline/20-mL HS vial
N8	0.8 μg gasoline/20-mL HS vial
N9	0.4 μg gasoline/20-mL HS vial

Sample collection and preparation (3/3)

- Simulated fire debris samples
 - Nylon carpet
 - Direct heat method




	Analyte concentration
FD1	100 μg gasoline + 0.25 g substrate/ 20-mL HS vial
FD2	50 μg gasoline+ 0.25 g substrate/ 20-mL HS vial
FD3	25 μg gasoline+ 0.25 g substrate/ 20-mL HS vial
FD4	12.5 μg gasoline+ 0.25 g substrate/ 20-mL HS vial
FD5	6.3 μg gasoline+ 0.25 g substrate/ 20-mL HS vial
FD6	3.1 μg gasoline+ 0.25 g substrate/ 20-mL HS vial
FD7	1.6 μg gasoline+ 0.25 g substrate/ 20-mL HS vial
FD8	0.8 μg gasoline+ 0.25 g substrate/ 20-mL HS vial
FD9	0.4 μg gasoline+ 0.25 g substrate/ 20-mL HS vial

Instrumental analysis

HS-SPME-GC/MS analysis




GC Oven Program Steps	Condition
GC oven initial temperature	40 °C
Hold time	5 min
Rate #1	10 °C/min
Oven temperature #1	150 °C
Rate #2	30 °C/min
Oven temperature #2	300 °C



HHS-SPME Steps	Condition		
Pre-Fiber Conditioning Temperature	250 °C		
Pre-Fiber Conditioning Time	60 s		
Pre-Incubation Time	300 s		
Incubation Temperature	80 °C		
Extraction Time	120 s		
Desorb to	GC injection port		
Desorption Time	120 s		
Post-Fiber Conditioning Temperature	250 °C		
Post-Fiber Conditioning Time	600 s		
GC Runtime	1200 s		

Image transformation

GC/MS data were transformed into 3 types of images



^{*}Note: All axes and labels of the images were removed during transfer learning.


Transfer learning (1/2)

Preparation of data sets

	Number of GC/MS data acquired				Total number of transformed images	
Training	Gasoline present	Neat samples	Brand A	63	315	39080% for training20% for validation
			Brand B	63		
			Brand C	63		
			Brand D	63		
			Brand E	63		
	Gasoline absent	Burned carpet		75	75	
Verification	Gasoline Neat sam		oles	90	180	195
	present	Simulated fire debris samples		90		
	Gasoline absent	Burned carpet		15	15	

Transfer learning (2/2)

Re-train GoogLeNet for gasoline detection





Predictions on the verification data set (1/5)

Training: validation accuracy

Extracted ion heatmap	Heatmap	TIC
100%	100%	100%


• Verification:






Predictions on the verification data set (2/5)

Comparison of prediction probability for simulated fire debris samples




Predictions on the verification data set (3/5)


Comparison of correct predictions for simulated fire debris samples



Predictions on the verification data set (4/5)


 Comparison of the predictions on the simulated fire debris samples at 1.6 - 100 μg gasoline sample/20-mL HS vial





Predictions on the verification data set (5/5)

 Comparison of classification performance between the extracted ion heatmap and four ML models



Conclusions

Experimental outcome

- TIC and heatmaps provided characteristic features of gasoline chemical profiles for transfer learning
- High performance for neat samples; limitation on fire debris samples
- Classification performance:
 - ➤ Heatmap > TIC
 - > Extracted ion heatmap > all ion range heatmap
 - Extracted ion heatmap > ML models

Intelligent workflow

- No dependency on manual feature extraction
- Achieved high accuracy without large-scale data collection
- More capable of discriminating mixtures compared to other ML models

Acknowledgements

This work was partly funded by the 2022-2023 Forensic Sciences Foundation (FSF) Lucas Research Grants. The opinions, findings, and conclusions or recommendations expressed in this presentation are those of the author(s) and do not necessarily reflect those of the FSF. The authors also appreciate the funding support from the Ministry of Education, Taiwan.



Thank you

Ting-Yu Huang txh038@shsu.edu