

Direct-to-PCR tissue preservation for DNA profiling

Amy Sorensen, MS¹; Clare Berry, BAS²; David Bruce, PhD³;
Michelle Gahan, PhD²; Sheree Hughes-Stamm, PhD¹; Dennis McNevin, PhD²

¹Department of Forensic Science, Sam Houston State University, USA

²National Centre for Forensic Studies, University of Canberra, Australia

³NSW Forensic and Analytical Science Service, Australia

68th Annual AAFS, Las Vegas, 2016



Sam Houston
State University



FSF Emerging Forensic Scientist Award
Oral Presentation

Warning!

There is some content in this presentation that some people may find disturbing.

DISCLOSURE

- Products used:
 - DNA Genotek Tissue Stabilizing Kit
 - DNAgard®
 - Quantifiler® Human DNA Quantification Kit
 - PowerPlex® 21
 - Globalfiler®

- The authors declare that they have no competing interests.



INTRODUCTION

- Disaster Victim Identification (DVI)
 - Mass fatalities
 - Limited/lack of facilities
 - Loss of electricity
 - Remote locations
 - Extreme temperatures & humidity
- Remote & rural forensic casework





INTRODUCTION



- Successful DNA typing is impacted by:
 - Adverse environmental conditions
 - Speed of collection
 - Preservation/storage of samples
 - Minimizes DNA degradation
 - Inhibits microbial activity
 - Room temp. storage
- Faster methods are desired
 - Tissue preservation releases DNA into solution ready for PCR



PREVIOUS LITERATURE

Allen-Hall, McNevin (2012)

- Fresh muscle
- Short-term storage (4-28 days)
 - 35°C, humidity 9-26%
- Results:
 - DNA recovered from preserved tissues
 - DNA in aliquots of preservatives
 - DESS, DNA Genotek, DNAgard, TENT
 - TENT failed to protect DNA from further degradation

Preservatives	Constituents
Dehydration	Oven drying at 35°C
Solid NaCl	Laboratory grade NaCl
DESS	20% DMSO, 0.25M EDTA, saturated with NaCl, pH 8.0
Ethanol	70% ethanol, 30% ddH ₂ O
Ethanol + EDTA	70% ethanol, 30% ddH ₂ O, 0.1mM EDTA
TENT buffer	10mM Tris, 10mM EDTA, 100mM NaCl, 2% Tween 20
RNAlater®	Proprietary
DNA Genotek Tissue Stabilizing Kit	Proprietary
DNAgard®	Proprietary



- Advantages to eliminating DNA extraction:

-


1. Preservatives contain high concentrations of PCR inhibitors
2. Components of tissue, by-products of decomposition, & contaminating material from soil also inhibit PCR



MATERIALS & METHODS

Preservative solutions that previously leached DNA into solution (2012)

- DMSO/EDTA saturated w/ salt(DESS)
- Tris, EDTA, NaCl, Tween 20 (TENT)
- DNAgard® [Biomatrica]
- DNA Genotek Tissue Stabilizing Kit [DNA Genotek]





MATERIALS & METHODS

Preservation of Tissue Samples

1. Fresh skeletal muscle (3 donors)

- 300mg preserved in 1 or 2mL of preservative
- Stored at 35°C for 3, 7, 14, & 28 days
- 20 or 50μL aliquot archived at -80°C for four years



UNIVERSITY OF
CANBERRA

2. Decomposing skin & muscle (2 cadavers)

- Tissue collected at 0, 6, 8, & 10 days of decomposition
- 30mg in 300μL of preservative
- Stored at 35°C for 1 month
- "No preservative" control



Sam Houston
State University



MATERIALS & METHODS


Days of Decomposition

Day 0

Day 6

Day 8

Day 10





MATERIALS & METHODS

Quantification of DNA in Preservatives:

- Quantifiler[®] Human DNA Quantification Kit [Thermo Fisher Scientific]
 - PCR inhibition measured by a delay in C_T for the IPC
 - When IPC was undetected → 1:10 or 1:20 dilution
 - Only DESS, DNAgard[®], & DNA Genotek from decomposed samples
- 7500 Real-Time PCR System [Thermo Fisher Scientific]



MATERIALS & METHODS

STR Genotyping

- Fresh tissue: 0.5ng of DNA added directly to PowerPlex® 21 System [Promega]
- Decomposed tissue: 0.5ng of DNA added directly to GlobalFiler® [Thermo Fisher Scientific]
- Samples amplified on 9700 & CE on 3500 [Thermo Fisher Scientific]

Statistical Analysis using SPSS [IBM]


- Shapiro-Wilk test for normality
- ANOVA (normally distributed data) or Kruskal-Wallis test (non-normal data)
- $p < 0.05$



RESULTS - FRESH

DNA Concentrations in Preservatives

- No PCR inhibition
- DNA yield highly dependent on donor tissue (Kruskal-Wallis p value=0.003)
 - Highest [DNA] from donor 1
- Significant difference between preservatives ($p=0.005$)
 - TENT yielding highest [DNA]
 - DNAgard® with the least
- [DNA] increases for longer storage (up to 28 days)






RESULTS - FRESH

STR Genotyping in Preservatives (PowerPlex 21[®])

- DESS & DNA Genotek generated full profiles
- DNAgard generated full profiles (except day 7)
- TENT produced profiles with 26-42 alleles, peak heights diminished at higher molecular weight loci (Penta D & E)
- Donor 2 had highest peak heights for all loci & all preservatives







RESULTS - DECOMPOSED

DNA Concentrations in Preservatives

Skin



Muscle






RESULTS - DECOMPOSED

DNA Concentrations in Preservatives

Skin




- Inhibition in preservatives surrounding decomp tissues
 - IPC only detected in TENT
 - 1:10 dil for DESS & DNAgard®
 - 1:20 dil for DNA Genotek
- All preservatives yielded more DNA than water control
- [DNA] increased from day 0 to 6, then decreased to nearly 0 at day 8 for B (skin & muscle) & at day 10 for A (consistent with time of bloat)




RESULTS - DECOMPOSED

STR Genotyping in Preservatives (GlobalFiler®)

Skin



Muscle






RESULTS - DECOMPOSED

STR Genotyping in Preservatives (GlobalFiler®)

Skin



- 1:10 dil for DESS & DNAgard, 1:20 dil for DNA Genotek required to relieve inhibition
- Full GF profiles from DESS, TENT, & DNA Genotek up to 8 days from cadaver A skin
- All preservatives yielded full profiles up to 6 days from cadaver A skin & cadaver B muscle



DISCUSSION

- Aliquots of preservative can be stored at -80°C for up to 4 years:
 - All except TENT are likely to produce full profiles
 - Long term archival does not diminish successful DNA typing
 - Storage at -80°C reduces PCR inhibition (when compared to samples stored at room temp.)
- DNA extraction can be eliminated
 - Faster DNA based DVI
 - Any increase in inhibition (from preservatives) can be easily diluted out
- Direct PCR may be suitable for DVI, provided that:
 - Tissue is collected before full bloat or entering active decay
 - DNA is of sufficient quantity & quality




CONCLUSIONS

- Decomposing tissue can be preserved at room temp.
- Tissue digestion & DNA extraction can be avoided
- Direct PCR approach for identifying fresh & decomposed tissue preserved at room temp. is possible
- Complete profiles can be obtained in a more timelier manner
- *Significant impact to address demands for immediate sample preservation & provide faster DNA identification*



FUTURE WORK

- Optimizing TENT buffer to better preserve DNA for storage & facilitate direct PCR
- Collection of DNA directly from decomposing bodies using FTA Elute cards, Whatman EasiCollect devices, Bode Buccal DNA Collection System, & traditional cotton swabs for room temp. storage



<http://www.mnready.com/buccal-006836%20card%20out.jpg>



<http://www.bodetech.com/wp-content/uploads/2011/01/bdcimg1.jpg>

ACKNOWLEDGEMENTS

- Southeast Texas Applied Forensic Science Facility (STAFS)
- Special thanks to the families of the ones that were donated
- Travel Funding provided by NIJ (GRF: 2015-R2-CX-0029)



Sam Houston
State University



REFERENCES

- Allen-Hall A, McNevin D. (2012) Human tissue preservation for disaster victim identification (DVI) in tropical climates. *Forensic Sci Int Genet* 6(5):653-657.

For more information, refer to:

- Sorensen A, Berry C, Bruce D, Gahan ME, Hughes-Stamm S, McNevin D. (2015) Direct-to-PCR tissue preservation for DNA profiling. *Int J Legal Med*



FSF Emerging Forensic Scientist Award
Oral Presentation

Questions?

Amy Sorensen
aes049@shsu.edu



Sam Houston
State University